

NMIN Capacity-Building Webinar

rhAmpSeq™ CRISPR: Multiplexed Amplicon Sequencing & Analysis



Adam Chernick

Field Application Manager
Integrated DNA Technologies (IDT)



rhAmpSeq™ CRISPR MULTIPLEXED AMPLICON SEQUENCING AND ANALYSIS

ADAM CHERNICK, PhD

November 16, 2020



The information in this presentation is confidential and proprietary. No part of the presentation may be used, reproduced or distributed without the express written permission of Integrated DNATechnologies, Inc.

THE WORLD'S LARGEST SUPPLIER OF CUSTOM NUCLEIC ACIDS



>64,000 oligos synthesized every day


- Major R&D teams (~100 employees)
- CRISPR
- NGS
- Genotyping
- qPCR
- Synthetic Biology
- RNAi
- Bioinformatics






IDT SUPPORTS THE CRISPR WORKFLOW

Alt-R CRISPR System – a complete workflow





CRISPR GENOME EDITING OVERVIEW





ON-TARGET AND OFF-TARGET EDITING

- Genome editing is not perfect
- CRISPR can result in edits throughout the genome
- On-target and off-target edits must be identified and quantified



WHAT DOES rhAmpSeq DO?

- Enables the best-in-class insertion/deletion (indel) quantification of genome editing
- Sequences on-target and off-target sites simultaneously
- Yields great uniformity and high “map rate” as well as quantitative readout
 - GUIDE-seq or other unbiased methods result in significant noise (>80% unmapped or non-target reads) in library prep; the rhAmpSeq system does not generate these wasted reads
- Provides detailed info at on-target sites, as well as off-target sites, all in one reaction
 - The condensed workflow saves time and reagents during design and experiment

The rhAmpSeq system is NOT an off-target identification method; it is a tool used to evaluate up to 5000 *known* or *predicted* off-target sites genome-wide.



STEP 1: IDENTIFY OFF-TARGET SITES

- Unbiased detection methods
 - Cell-based—GUIDE-seq
 - *In vitro*—SITE-Seq
- *In silico* prediction methods
 - IDT “gRNA checker”
 - Cosmid, CRISTA, RGEN

GUIDE-seq results

Internal, improved protocols to share with customers
 Tsai, *et al.*, *Nat Biotech* (2015)

Mismatch	Assay	Chr	Guide Coordinates		Strand	Sequence Alignments	GS Reads
			Start	Stop		G T T G G A G C A T C T G A G T C C A G	
0	AR_iGS_1	chrX	67545904	67545924	+	G T T G G A G C A T C T G A G T C C A G	4030
4	AR_iGS_2	chr17	14626784	14626804	-	A A T G G G G G C A T C T G A G T C C A T	1671
3	AR_iMS_2	chr1	27592693	27592713	-	G A T G G A G C A A C C G A G T C C A G	1190
2	AR_iGS_3	chr7	22126335	22126355	-	A T T G G A G C C T C T G A G T C C A G	940
2	AR_iGS_4	chr20	46362541	46362561	-	G T T G G A G A A A C T G A G T C C A G	562
1	AR_iMS_3	chr8	70014903	70014923	-	G T T G G A G C C T C T G A G T C C A G	418
3	AR_iMS_4	chr12	122113358	122113378	-	G G A G G A G C A C C T G A G T C C A G	409
3	AR_iGS_5	chr10	75605948	75605968	-	G G A A G A G C A T C T G A G T C C A G	341

STEP 2: STUDY THEM IN DETAIL WITH rhAmpSeq




STEP 1: IDENTIFICATION

- Unbiased detection methods
 - Cell-based—GUIDE-seq
 - *In vitro*—SITE-Seq

GUIDE-seq results

Internal, improved protocols to share with customers
 Tsai, *et al.*, *Nat Biotech* (2014)



TESTS

prediction methods
 “gRNA checker”
 CRISPRid, CRISTA, RGEN

Sequence Alignments		GS Reads
end	G T T G G A G C A T C T G A G T C C A G	4030
	G T T G G A G C A T C T G A G T C C A G	1671
	A A T G G G G G C A T C T G A G T C C A T	1190
	G A T G G A G C A A C C G A G T C C A G	940
	A T T G G A G C C T C T G A G T C C A G	562
	G T T G G A G A A A C T G A G T C C A G	418
	G G A G G A G C A C C T G A G T C C A G	409
	G G A A G A G C A T C T G A G T C C A G	341

STEP 2: STUDY


VITH rhAmpSeq




rhAmpSeq METHODOLOGY?

Step 1: rhPCR (1.2 hrs)

Activation of rhPrimers by
RNase H2 Cleavage




Amplification by DNA
Polymerase



Step 2: Universal PCR (1 hr)


Index rhPCR amplicons
with P7 & P5 universal
primers



Quantify each rhAmpSeq library
for normalization and sample pooling



rhAmpSeq LIBRARY FORMATION



LSS: Locus-specific sequence

Assay: A pair of forward and reverse rhAmp primers designed to amplify a specific region

Panel: Collection of assays (forward panel and reverse panel)



rhAmp PCR TECHNOLOGY INNOVATED AT IDT

BMC Biotechnology (2011) 11:80

METHODOLOGY ARTICLE

Open Access

RNase H-dependent PCR (rhPCR): improved specificity and single nucleotide polymorphism detection using blocked cleavable primers


Joseph R Dobosy, Scott D Rose, Kristin R Beltz, Susan M Rupp, Kristy M Powers, Mark A Behlke* and Joseph A Walder

Blocked inactive primers anneal to target

RNase H2 recognizes hybridized internal RNA base

RNase H2 cleaves hybridized primers

DNA polymerase extends newly unblocked primers






BENEFITS OF rhAmp PCR


RNase H2 is sensitive to mismatches around the RNA base cleavage site, thus enabling:

- Higher specificity
- Lower potential for primer-dimers
- Higher multiplexing capability
- Modular panel designs
- Flexible amplicon sizes





rhAmpSeq DESIGN APPROACH




Design

1. Flexible amplicon size
2. Exhaustive primer design around each target site
3. Comprehensive QC



rhAmpSeq DESIGN APPROACH



PUBLICATIONS: ON/OFF-TARGET WORKFLOW



Molecular Therapy
Methods & Clinical Development
Original Article



Increasing CRISPR Efficiency and Measuring Its Specificity in HSPCs Using a Clinically Relevant System

Jenny Shapiro,^{1,5} Ortal Iancu,^{1,5} Ashley M. Jacobi,² Matthew S. McNeill,² Rolf Turk,² Garrett R. Rettig,² Ido Amit,³ Adi Tovin-Recht,¹ Zohar Yakhini,^{3,4} Mark A. Behlke,² and Ayal Hendel¹

¹Institute of Nanotechnology and Advanced Materials, The Mina and Everard Goodman Faculty of Life Sciences, Bar-Ilan University, Ramat-Gan 5290002, Israel; ²Integrated DNA Technologies, Coralville, IA 52241, USA; ³Department of Computer Science, Interdisciplinary Center, Herzliya 4610101, Israel; ⁴Department of Computer Science, Technion-Israel Institute of Technology, Haifa 3200003, Israel



bioRxiv

THE PREPRINT SERVER FOR BIOLOGY

CIS checkpoint deletion enhances the fitness of cord blood derived natural killer cells transduced with a chimeric antigen receptor

May Daher, Rafet Basar, Elif Gokdemir, Natalia Baran, Nadima Uprety, Ana Karen Nunez Cortes, Mayela Mendt, Lucila Nassif Kerbauy, Pinaki P. Banerjee, Mayra Hernandez Sanabria, Nobuhiko Imahashi, Li Li, Francesca Lorraine Wei Inng Lim, Mohsen Fathi, Ali Rezvan, Vakul Mohanty, Yifei Shen, Hila Shaim, Junjun Lu, Gonca Ozcan, Emily Ensley, Mecit Kaplan, Vandana Nandivada, Mustafa Bdaawi, Sunil Acharya, Yuanxin Xi, Xinhai Wan, Duncan Mak, Enli Liu, Sonny Ang, Luis Muniz-Feliciano, Ye Li, Jing Wang, Shahram Kordasti, Nedyalko Petrov, Navin Varadarajan, David Marin, Lorenzo Brunetti, Richard J. Skinner, Shangrong Lyu, Leiser Silva, Rolf Turk, Mollie S. Schubert, Garrett R. Rettig, Matthew S. McNeill, Gavin Kurgan, Mark A. Behlke, Heng Li, Natalie W. Fowlkes, Ken Chen, Marina Konopleva, Richard Champlin, Elizabeth J. Shpall, Katayoun Rezvani

doi: <https://doi.org/10.1101/2020.03.29.014472>

This article is a preprint and has not been certified by peer review [what does this mean?].

REGULAR ARTICLE

blood advances

Large-scale GMP-compliant CRISPR-Cas9-mediated deletion of the glucocorticoid receptor in multivirus-specific T cells

Rafet Basar,^{1,*} May Daher,^{1,*} Nadima Uprety,^{1,*} Elif Gokdemir,¹ Abdullah Alsuliman,¹ Emily Ensley,¹ Gonca Ozcan,¹ Mayela Mendt,¹ Mayra Hernandez Sanabria,¹ Lucila Nassif Kerbauy,^{2,3} Ana Karen Nunez Cortes,¹ Li Li,¹ Pinaki P. Banerjee,¹ Luis Muniz-Feliciano,¹ Sunil Acharya,¹ Natalie W. Fowlkes,⁴ Junjun Lu,¹ Sufang Li,¹ Stephan Mielke,⁵ Mecit Kaplan,¹ Vandana Nandivada,¹ Mustafa Bdaawi,¹ Alexander D. Kontoyiannis,⁶ Ye Li,¹ Enli Liu,¹ Sonny Ang,¹ David Marin,¹ Lorenzo Brunetti,^{7,8} Michael C. Gundry,^{7,8} Rolf Turk,⁹ Mollie S. Schubert,⁹ Garrett R. Rettig,⁹ Matthew S. McNeill,⁹ Gavin Kurgan,⁹ Mark A. Behlke,⁹ Richard Champlin,¹ Elizabeth J. Shpall,¹ and Katayoun Rezvani¹

Published online 31 May 2019

Nucleic Acids Research, 2019, Vol. 47, No. 15 7955–7972
doi: 10.1093/nar/gkz475

Highly efficient editing of the β -globin gene in patient-derived hematopoietic stem and progenitor cells to treat sickle cell disease

So Hyun Park^{1,†}, Ciaran M. Lee^{1,†}, Daniel P. Dever², Timothy H. Davis¹, Joab Camarena², Waracharee Srifa², Yankai Zhang³, Alireza Paikari³, Alicia K. Chang³, Matthew H. Porteus², Vivien A. Sheehan^{3,*} and Gang Bao^{1,*}

nature
medicine

ARTICLES

<https://doi.org/10.1038/s41591-018-0137-0>

A high-fidelity Cas9 mutant delivered as a ribonucleoprotein complex enables efficient gene editing in human hematopoietic stem and progenitor cells

Christopher A. Vakulskas^{1,7}, Daniel P. Dever^{2,7}, Garrett R. Rettig¹, Rolf Turk¹, Ashley M. Jacobi¹, Michael A. Collingwood¹, Nicole M. Bode¹, Matthew S. McNeill¹, Shuqi Yan¹, Joab Camarena², Ciaran M. Lee³, So Hyun Park³, Volker Wiebking², Rasmus O. Bak^{4,5}, Natalia Gomez-Ospina², Mara Pavel-Dinu², Wenchen Sun⁶, Gang Bao³, Matthew H. Porteus² and Mark A. Behlke¹

nature
biotechnology

ARTICLES

<https://doi.org/10.1038/s41587-020-0555-7>



CHANGE-seq reveals genetic and epigenetic effects on CRISPR-Cas9 genome-wide activity

Cicera R. Lazzarotto¹, Nikolay L. Malinin¹, Yichao Li¹, Ruochi Zhang², Yang Yang², GaHyun Lee¹, Eleanor Cowley³, Yanghua He^{1,8}, Xin Lan¹, Kasey Jividen¹, Varun Katta¹, Natalia G. Kolmakova⁴, Christopher T. Petersen⁵, Qian Qi¹, Evgheni Strelcov^{6,7}, Samantha Maragh⁴, Giedre Krenciute⁵, Jian Ma², Yong Cheng¹ and Shengdar Q. Tsai^{1,8}

IDT
INTEGRATED DNA TECHNOLOGIES

PUBLICATIONS: ON/OFF-TARGET WORKFLOW



Molecular Therapy
Methods & Clinical Development
Original Article



Increasing CRISPR Efficiency and Measuring Its Specificity in HSPCs Using a Clinically Relevant System

Jenny Shapiro,^{1,5} Ortal Iancu,^{1,5} Ashley M. Jacobi,² Matthew S. McNeill,² Rolf Turk,² Garrett R. Rettig,² Ido Amit,³ Adi Tovin-Recht,¹ Zohar Yakhini,^{3,4} Mark A. Behlke,² and Ayal Hendel¹

¹Institute of Nanotechnology and Advanced Materials, The Mina and Everard Goodman Faculty of Life Sciences, Bar-Ilan University, Ramat-Gan 5290002, Israel; ²Integrated DNA Technologies, Corvallis, IA 52241, USA; ³Department of Computer Science, Interdisciplinary Center, Herzliya 4610101, Israel; ⁴Department of Computer Science, Technion-Israel Institute of Technology, Haifa 3200003, Israel



bioRxiv
THE PREPRINT SERVER FOR BIOLOGY

CIS checkpoint deletion enhances the fitness of cord blood derived natural killer cells transduced with a chimeric antigen receptor

May Daher, Rafet Basar, Elif Gokdemir, Natalia Baran, Nadima Uprety, Ana Karen Nunez Cortes, Mayela Mend Lucila Nassif Kerbauy, Pinaki P. Banerjee, Mayra Hernandez Sanabria, Nobuhiko Imahashi, Li Li, Francesca Lorraine Wei Inng Lim, Mohsen Fathi, Ali Rezvan, Vakul Mohanty, Yifei Shen, Hila Shaim, Junjun Lu, Gonca Ozcan, Emily Ensley, Mecit Kaplan, Vandana Nandivada, Mustafa Bdaawi, Sunil Acharya, Yuanxin Xi, Xinhai Wan, Duncan Mak, Enli Liu, Sonny Ang, Luis Muniz-Feliciano, Ye Li, Jing Wang, Shahram Kordasti, Nedyalko Petrov, Navin Varadarajan, David Marin, Lorenzo Brunetti, Richard J. Skinner, Shangrong Lyu, Leiser Silva, Rolf Turk, Mollie S. Schubert, Garrett R. Rettig, Matthew S. McNeill, Gavin Kurgan, Mark A. Behlke, Heng Li, Natalie W. Fowlkes, Ken Chen, Marina Konopleva, Richard Champlin, Elizabeth J. Shpall, Katayoun Rezvani

doi: <https://doi.org/10.1101/2020.03.29.014472>

This article is a preprint and has not been certified by peer review [what does this mean?].

REGULAR ARTICLE

blood advances

Large-scale GMP-compliant CRISPR-Cas9-mediated deletion of the glucocorticoid receptor in multivirus-specific T cells

Rafet Basar,^{1,*} May Daher,^{1,*} Nadima Uprety,^{1,*} Elif Gokdemir,¹ Abdullah Alsuliman,¹ Emily Ensley,¹ Gonca Ozcan,¹ Mayela Mendt,¹ Mayra Hernandez Sanabria,¹ Lucila Nassif Kerbauy,^{2,3} Ana Karen Nunez Cortes,¹ Li Li,¹ Pinaki P. Banerjee,¹ Luis Muniz-Feliciano,¹ Sunil Acharya,¹ Natalie W. Fowlkes,⁴ Junjun Lu,¹ Sufang Li,¹ Stephan Mielke,⁵ Mecit Kaplan,¹ Vandana Nandivada,¹ Mustafa Bdaawi,¹ Alexander D. Kontoyiannis,⁶ Ye Li,¹ Enli Liu,¹ Sonny Ang,¹ David Marin,¹ Lorenzo Brunetti,^{7,8} Michael C. Gundry,^{7,8} Rolf Turk,⁹ Mollie S. Schubert,⁹ Garrett R. Rettig,⁹ Matthew S. McNeill,⁹ Gavin Kurgan,⁹ Mark A. Behlke,⁹ Richard Champlin,¹ Elizabeth J. Shpall,¹ and Katayoun Rezvani¹

Published online 31 May 2019

Nucleic Acids Research, 2019, Vol. 47, No. 15 7955–7972
doi: 10.1093/nar/gkz475

Highly efficient editing of the β -globin gene in patient-derived hematopoietic stem and progenitor cells to treat sickle cell disease

So Hyun Park^{1,†}, Ciaran M. Lee^{1,†}, Daniel P. Dever², Timothy H. Davis¹, Joab Camarena², Waracharee Srifa², Yankai Zhang², Alireza Paikari³, Alicia K. Chang³, Matthew H. Porteus², Vivien A. Sheehan² and Gang Bao^{1,*}

nature
medicine

ARTICLES

<https://doi.org/10.1038/s41591-018-0137-0>

A high-fidelity Cas9 mutant delivered as a ribonucleoprotein complex enables efficient gene editing in human hematopoietic stem and progenitor cells

Christopher A. Vakulskas^{1,7}, Daniel P. Dever^{2,7}, Garrett R. Rettig¹, Rolf Turk¹, Ashley M. Jacobi¹, Michael A. Collingwood¹, Nicole M. Bode¹, Matthew S. McNeill¹, Shuqi Yan¹, Joab Camarena², Ciaran M. Lee³, So Hyun Park³, Volker Wiebking², Rasmus O. Bak^{4,5}, Natalia Gomez-Ospina², Mara Pavel-Dinu², Wenchen Sun⁶, Gang Bao³, Matthew H. Porteus² and Mark A. Behlke¹

nature
biotechnology

ARTICLES

<https://doi.org/10.1038/s41587-020-0555-7>




CHANGE-seq reveals genetic and epigenetic effects on CRISPR-Cas9 genome-wide activity

Cicera R. Lazzarotto¹, Nikolay L. Malinin¹, Yichao Li¹, Ruochi Zhang², Yang Yang², GaHyun Lee¹, Eleanor Cowley³, Yanghua He^{1,8}, Xin Lan¹, Kasey Jividen¹, Varun Katta¹, Natalia G. Kolmakova⁴, Christopher T. Petersen⁵, Qian Qi¹, Evgheni Strelcov^{6,7}, Samantha Maragh⁴, Giedre Krenciute⁵, Jian Ma², Yong Cheng¹ and Shengdar Q. Tsai^{1,8}

IDT
INTEGRATED DNA TECHNOLOGIES



rhAmpSeq TECHNOLOGY AND Alt-R HiFi Cas9



rhAmpSeq system savings:

- 8 x 40 assay designs
- ~1500 individual PCRs reduced to <96
 - Master mix
 - gDNA
- Library quantification
- Full-time equivalent hours
 - Months to days




DATA ANALYSIS PIPELINE

Off-target identification
(stable Cas9, modified gRNAs)




rhAmpSeq design,
amplification,
and sequencing
(multiple conditions,
high-throughput)



The screenshot displays the IDT web interface on the left and a Windows File Explorer window on the right. The IDT interface includes a navigation bar with 'Files', 'Samples', 'Analyses', 'Aggregations', and 'Queries'. Below this is a 'Files' section with a 'Choose a file or drag it here' area and a table with columns for 'File Name', 'Creation Date', 'Size', and 'Source'. The Windows File Explorer window shows a folder named 'alpha-test' containing several files, including 'Singleplex_GeneX_Control_R1.fastq.gz', 'Singleplex_GeneX_Control_R2.fastq.gz', 'Singleplex_GeneX_EditedHDR_R1.fastq.gz', 'Singleplex_GeneX_EditedHDR_R2.fastq.gz', 'Singleplex_GeneX_EditedOnly_R1.fastq.gz', and 'Singleplex_GeneX_EditedOnly_R2.fastq.gz'. A red circle highlights these files, and a red arrow points from the 'Choose a file or drag it here' area in the IDT interface to the highlighted files in the File Explorer.





CRISPAItRations FOR THE ANALYSIS OF rhAmpSeq DATA





COMPARATIVE ANALYSIS OF CRISPR-NGS TOOLS





CONCLUSIONS

- rhAmp PCR enables high multiplexing levels.
- Up to 5000 targets are designed for a single-pool reaction.
- Simultaneously quantify editing for on- and off-target loci for CRISPR experiments.
- Streamlined, simple, 2-step PCR protocol generates NGS-ready libraries in approximately 4 hours.
- High on-target rate and coverage uniformity from as little as 10 ng of DNA input.
 - Increased % on-target reads
 - Higher degree of multiplexing
 - Better uniformity in coverage across the panel



THANK YOU

www.idtdna.com/CRISPR
CRISPR@idtdna.com
achernick@idtdna.com



NMIN Capacity-Building Webinar

Q&A



Adam Chernick

Field Application Manager
Integrated DNA Technologies (IDT)



Nanomedicines Innovation Network (NMIN)



<https://nanomedicines.ca>



info@nanomedicines.ca



[@NMIN_NCE](https://twitter.com/NMIN_NCE)

